



5905 Smiley Drive
 Culver City, CA 90232
 PH: 310-841-0123
 FAX: 310-841-0125

BACKING HANDLING AND DISINFECTING PROCEDURES

JC Backings employees are required to wear face masks while on site, socially distance throughout the day and wash their hands or use hand sanitizer every hour and after handling any items/backing(s).

| Type of Backing | JC Backings Employee Protocol |
|--|--|
| Photo and Digital Print Backing Returns Within LA County | <ul style="list-style-type: none"> • Off-load backing(s) in the parking lot with forklift or use hands with horns when necessary. • Wash or sanitize your hands before spraying tube(s). • Spray the tube with CDC-approved disinfectant. • Bring into warehouse and quarantine for 48 hours. • Backing(s) will be checked after 48 hours. |
| Photo and Digital Print Backing Returns from Outside of LA County and Out-of-State | <ul style="list-style-type: none"> • Off-load backing(s) in the parking lot with forklift or use hands with horns when necessary. • Wash or sanitize your hands before spraying tube(s). • Spray the plastic wrapping with CDC-approved disinfectant. • Carefully cut the plastic wrapping off. • Spray the tube with CDC-approved disinfectant. • Bring into warehouse and quarantine for 48 hours. • Backing(s) will be checked after 48 hours. |
| Painted Backing Returns Within LA County | <ul style="list-style-type: none"> • Off-load backing(s) in the parking lot with forklift or use hands with horns when necessary. • Wash or sanitize your hands before wrapping backing(s) in plastic. • Spray the plastic wrapped backing(s) with CDC-approved disinfectant. • Once dry bring plastic wrapped backing(s) into warehouse and quarantine for 48 hours. • Backing(s) will be checked after 48 hours. |
| Painted Backing Returns from Outside of LA County and Out-of-State | <ul style="list-style-type: none"> • Off-load backing(s) in the parking lot with forklift or use hands with horns when necessary. • Wash or sanitize your hands before spraying tube(s)/crate(s). • Spray the plastic wrapping (if tubed) or the crate with CDC-approved disinfectant. • Carefully cut the plastic wrapping off the tube. • Spray the tube with CDC-approved disinfectant. • Bring tube and/or crate into warehouse and quarantine for 48 hours. • Backing(s) will be checked after 48 hours. |

**Thank you for your patience and cooperation with our requirements.
 If you have questions, please call our office at 310-841-0123.**